

Everyday care can be wrapped
with smile



 wrappon

Auto-wrapping (tear-and seal) toilet

 wrappon

 wrappon · AILE  wrappon · EVERI



■ Contact Information

 NIHON SAFETY CO.,LTD.

<http://www.wrappon.com>


Wrappon Division, Nihon Safety Co., Ltd.

11F, Ichibancho Tokyu Bldg.

21 Ichibancho, Chiyoda-ku, Tokyo 102-0082 JAPAN

Phone : +81-3-6369-2223

E-mail: kaigo@nihonsafety.com

 wrappon is a registered trademark of NIHONSAFETY CO.,LTD.

These product descriptions are effective as of September 2016. For product improvement purposes, the specifications and external appearance may be changed without prior notice.

■ Agent



佳澤工程有限公司
Delta Pyramax Engineering Ltd.

3/F, Unit A, Kader Building, 22 Kai Cheung Road, Kowloon Bay, HK
Phone +852 2511 2118 | joyceli@dpx.hk | www.deltapyramax.com

WPAEn1609-02

 NIHON SAFETY CO.,LTD.

Basic Type wrappon · EVERI

4th Eldercare
Innovation
Awards 2016



Spacious wooden chair type can also be used for well-built persons. This auto-wrapping model is a basic type with fundamental functions such as adjustable heights and detachable armrests.

Adjustable armrest for easier movement.
For easier movement from bed to the toilet, left and right armrests can be detachable.



Height adjustable Armrests.
Two heights can be set per your needs.



Soft toilet seat.
For persons who sit on longer time and have smaller hips.



Three height adjustable legs.
To meet with the length and preference, legs height can be adjusted three levels with 3cm increment.



How to use the auto-wrapping defecation treatment toilet.

It is the same way for both Everi and Aile. Just push the flush buttons so anybody can use this toilet.



Pour the designated coagulant.
To coagulate the moisture out of waste.



Push the flush button.
After relieving yourself, auto-wrapping starts.



Pick up the wrapped film
Cut off is completed when Beep Sounds heard. Then, pick it up.

**Beep!
Beep!
Beep!**

Smart Type *New Release* wrappon · AILE

wrappon × Richell

To meet with customer needs, we improved the wooden chair. This is a new collaboration between compact size Richell wooden chair which is powered by Wrappon auto-wrapping defecation treatment unit.

Three different colors



Beige



Navy



Yellow

Flippable armrests.
Flippable armrest is easier to move from the bed to the toilet horizontally..



Three height adjustable armrests.
Three heights can be adjusted per your size and conditions(18, 21, 24cm).



Easy movement.
Casters installed for easier movement.



Easy access from the back.
Toilet seat cover can be falled down backward so easy access from the back.




Toilet seat height adjustment.
Three different heights can be set.(40, 43, 46cm)



Change to Bucket type.
Easy change from auto-wrapping to the bucket cleaning.



Selectable toilet seat.
Can change toilet seat to soft toilet seat and to heated toilet seat as option.



Peace of mind when blackout.
You can use this toilet when there is blackout as you can use rechargeable battery and DC cable as option.



In an elderly care life, great distresses of defecation are smell, poor hygiene and a chore of bucket cleaning. Burdens are increasing more and more for not only a person himself/herself, but also his/her family as well. Wrappon can contribute to reduce the burdens for both care receiver and care giver and create mentally healthy environment.

In 2015, Wrappon was selected as Robot care development and introduction project

Things that you can change by Wrappon

Bucket cleaning



Bucket cleaning is not required as auto-wrapping defecation treatment unit is installed. It wraps human waste automatically. The wrapped film is sealed and cut and drop the wrapped film on tray.

Ageing-in-place



Simple operation. Pour the coagulant and push the button. The rest is done by the unit so this helps the care receiver to live independently.

Running cost



Consumables cost is 50 yen per toilet. Normal portable toilet requires bucket cleaning with coagulant and 60 L water. Plus the burden of bucket cleaning chores.

No smell



By sealing the film, the smell contains inside and prevent the smell going outside the room.

Reduce the burden



No worry his/her waste to be seen by other persons. No worry to ask for cleaning. Defecation care can be done by himself/herself with pride.

Family unity




Family unity can be broken by toilet cleaning. We have received from users who regained the family unity by using Wrappon.

What is auto-wrapping defecation treatment unit?

A bucket cleaning is not needed as wrapping waste and smell by heat sealed automatically.

Functions of auto-wrapping

1. Starting to wrap



After defecation, start auto-wrapping by just pushing the button.

2. Film feeding



Enclosed not only waste but also soiled area into the wrapped film.

3. Thermo-bonding



Thermo-bonding to conceal waste.

4. Cut off



After two minutes, the film is cut off automatically.

Easy to operate


After two minutes, the film is cut.

Thermo-bonding

Smell is also treated

Enclosed not only waste but also soiled area into the wrapped film.

Beautifully treated



Voices from end-users

We have received many happy voices from our end-users. Smiles are spreading nationwide by Wrappon users.

Without Wrappon Everi, I get to used to the smelling environment as normal condition.

Comfort care by Wrappon

Hyogo Prefecture Ms. Yamamoto The user: her mother(80's)

"Wrappon is easy to use as no smell and easy treatment after defecation. If I used the other portable toilet, I had a big burden not only for the smell but also for my care life. Since I work daytime, I cannot take care of the bucket cleaning in a daytime. It would be a big burden for my mother to keep the human waste until I go home to get rid of these hygienically. She keeps her clothes in the room. It is not recommended for her to visit Day care service wearing those smelly clothes. I would be very sad if she went to visi Day care service wearing smelly clothes. Without Wrappon, it would be a standard life for my mother with those smelly clothes."



It is much easier to get rid of human waste with Wrappon and it is changed my image of defecation care.

Comfort care by Wrappon



Shizuoka Prefecture Ms. Suzuki The user: his mother(70's)

"My mother was used to take care of her mother. At that time, the bucket cleaning toilet was used so toilet has always been difficult to deal with as daily chores. But, after using Wrappon, toilet has been changed to be happy and easy. This is totally different image that I used to have.

My mother is very satisfied as no smell spreading into the room. I'm a little concerned with the defecation care, but it is cleared that it can be treated same as normal garbage



SPECIFICATIONS

PRODUCT

wrappon • Eviri

Product Number	EVSE0003JH	JAN Code	4589922490015
TAIS Code	00942-0000009		
Sizes	Main body : W 565mm, D 560mm, H 740・770・800mm Toilet seat height: 370mm, 400mm, 430mm Armrest height from the toilet seat : 210mm, 240mm Seat width : 470mm		
Weight	Main body approximately 21kg		
Materials	Chair : Rubber wood with urethane coating Toilet seat : Polypropylene, Polyurethane		
Power supply	AC100V±10V	Power consumption	75w

※ラップボンの使用には、必ず消耗品が必要となります。



4th Eldercare Innovation Awards 2016

wrappon • Aile New Release

Product Number(Beige)	ALSE002BJH	JAN Code(Beige)	4589922490145
Product Number(Navy)	ALSE002NJH	JAN Code(Navy)	4589922490152
Product Number(Yellow)	ALSE002YJH	JAN Code(Yellow)	4589922490169
TAIS Code	00942-000011		
Sizes	Main body : W 508mm, D 530mm, H 830・860・890mm Toilet seat height: 400mm, 430mm, 460mm Armrest height from the toilet seat : 180,210, 240mm Seat width : 400mm		
Weight	Main body approximately 25kg		
Materials	Chair : Rubber wood with urethane coating Toilet seat : Polypropylene		
Power supply	AC100V±10V	Power consumption	75w

※ラップボンの使用には、必ず消耗品が必要となります。



Eviri ...For Wrappon・Everi Aile ...For Wrappon・Aile

Consumables

WRAPPON Designated Consumable Set for elder

Eviri Aile

Product Number	EVSEW502J
JAN Code	4589922490022
Content of set	Designated paper coagulant, film cassette type 2

*You need to place the order for the first time with the product.

Designated paper coagulant

Product Number	COCOSTP1J
JAN Code	4589922490046
Material	Recycled paper
Capacity	6L(60 times)

Film cassette Type2

Product Number	C0F1C260J
JAN Code	4589922490039
Material	Polyethylene film
Capacity	Approximately 60 times



* You need to order the above when using.

Option

Soft toilet seat (Replacement) Aile

Product Number	RCTSS51J
JAN Code	4589922490176
Material	Body : EVA resin Core plate : Polypropylene



Heated seat(replacement) Aile

Product Number	RCT5WS01J
JAN Code	4589922490183
Material	ABS resin
Safety device	Electrically controlled・Temperature Fuse(98℃)
Rated power	AC100V
Power consumption	42 W
Power cord	Main body : 1.2m Power source : 1.6m
Surface temperature	Controlled temperature at 20℃ from 20℃~40℃



Toilet paper holder(replacement) Aile

Product Number	ALPHR0015
JAN Code	4589922490190
Material	Toilet paper holder (partial) : Poly Paper holder : polypropylene Fixing jig : Rubber wood



車用DCケーブル Bタイプ Aile

品番	COCCB10MS
JANコード	4589922490268
寸法	10m
規格	12V車用 電源スイッチ付き



Designated Battery Aile

Product Number	BALPL5DXJH
JAN Code	4589922490251
Material	ABS resin
Battery capacity	DCSV Max 2A
USB output	DC12V
Output voltage	DC12V
Rechargeable battery	Lithium ion
Charging time	Approximately 4 hours 70 times (fully charged)

